PCD224-C

Process Challenge Device

Test Pack for Steam sterilization processes.

Usage

For routine load release, especially for loads containing implants. For routine sterilizer monitoring. For periodic validation of the sterilizer (after installation, repair, relocation).

Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity displacement steam autoclaves between 132-135 °.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11140-1:2014, ISO 11138-1:2017, ISO 11138-3:2017, ANSI/AAMI ST79.

Classification

Class 1, according to risk.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

FDA 510(k)

K191021

Characteristics

Porous cards system holding a Bionova® BT224 Self-contained Biological Indicator and a moving-front IT26-C Chemical Integrator (Type 5 Integrating Indicator according to ISO 11140-1: 2014). The system also contains a self-adhesive Record Card, PCDBI-C-RC, where data on sterilization cycle may be written.

Cardboard box (10.5 cm x 6.5 cm x 1.8 cm approximately) with a Process Indicator printed (Type 1 according to ISO 11140-1: 2014).

Process Indicator:

Initial color: **blue**

Final color: **dark grey or black**

IT26-C Chemical Integrator (within the challenge device):

The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met.

Integration condition is calibrated against the kill time of a 10^6 *G. stearothermophilus* ATCC 7953 spores' population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

BT224 Ultra Rapid Self-contained Biological Indicator (within the challenge device): $\geq 10^6$ *G. stearothermophilus* ATCC 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

Final results: fluorescence readout after 20-minutes incubation at 60 °C (sensitivity \geq 97%). A visual confirmation of the result is optional (color change of the medium caused by pH change) and can be made after incubation at 60 °C for 48 hours and/or for 7 days.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

Environmental conditions during manufacture

T= 15-30 °C, RH= 30-80 %.

Storage conditions

T= 10-30 °C, RH= 30-80 %, keep out of direct light.

Transport conditions

Storage conditions should be strictly followed. Carry in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf life

2 years.

Packing

Packing information: product code and description, regulation, process for intended use, conditions for use, instructions for use, storage conditions, manufacturer information and data on box's label.

Labelling

On product's packing: product code, batch number, batch number of the Self-Contained Biological Indicator, manufacture date, expiration date, regulation, bar code and datamatrix code.

On the BT224 Self-Contained Biological Indicator: label printed with Steam reactive ink (turns to brown), product's code, process for intended use, name of organism, batch number and expiration date.

On the Record Card, PCDBI-C-RC: product's code, description and process for intended use

On the moving-front Chemical Integrator IT26-C: product's code, process for intended use and classification according to regulation.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

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