Infection Control Division | **Process Challenge Devices**

KPCD224-C Process Challenge Device Kit

For steam sterilization processes

Usage

United States

Terragene® Bionova® PCD224-C provides a defined challenge resistance against the claimed cycles shown below and demonstrated resistance equivalence to the ANSI/AAMI 16 towel pack. The device provides routine monitoring and sterilizer qualification testing steam sterilization processes.

Autoclave/Steam cycles	
Gravity displacement	
132 °C	135 °C
10 min, 15 min, 25 min	10 min
Dynamic air removal (vacuum assist)	
132 °C	135 °C
4 min	3 min
Fluorescence read time	pH color change
20 min	48 hours

Outside the United States

Terragene® Bionova® PCD224-C has been designed for a rapid and easy monitoring of steam sterilization processes at 132-135 °C ≥ 4 minutes.

Applicable regulation

Designed under Quality Management Systems standards ISO 13485:2016 / EN ISO 13485:2016.

ISO 11138-1:2017, ISO 11138-3:2017, ISO 11138-8:2021, ISO 11140-1:2014, ANSI/AAMI ST79:2017.

Authorization

ANMAT (Argentinian National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Classification

Class 1, according to risk (ANMAT).

FDA 510(k)

K191021.

Characteristics

The kit consists of:

- 25 PCD224-C: Each of these porous card systems consists of a cardboard box (10.6 cm x 6.8 cm x 1.7 cm) with a Type 1 process indicator, containing:
 - A Self-Contained Biological Indicator Bionova® BT224.

Infection Control Division | Process Challenge Devices

KPCD224-C Process Challenge Device Kit

For steam sterilization processes

- A Type 5 chemical indicator Integron® IT26-C (moving-front integrator).

- A self-adhesive record card, PCDBI-C-RC, where data on sterilization cycle may be written.

- 25 Bionova® BT224 Biological Indicators (to be used as a positive controls).

PCD Process indicator:

Initial color: **blue**.

Final color: **dark grey or black**.

IT26-C integrator indicator (within the PCD):

The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met.

Integration condition is calibrated against the kill time of a 10^6 *Geobacillus stearothermophilus* ATCC®* 7953 spore population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: Saturated steam at 121 °C, 128 °C, 135 °C.

BT224 Self-contained Ultra Rapid Biological Indicator (within the PCD):

$\geq 10^6$ *G. stearothermophilus* ATCC® 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with steam reactive ink (color change: pink to brown).

Final results: By fluorescence after 20-minute incubation at 60 °C (sensitivity $\geq 97\%$). An optional visual pH color change confirmation could be made after 48 hours of incubation. A 7-day readout is optional to confirm the 20-minute result by color change.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

*ATCC® is a registered trademark of American Type Culture Collection.

Infection Control Division | **Process Challenge Devices**

KPCD224-C Process Challenge Device Kit

For steam sterilization processes

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep in a dark place.

Transport conditions

It is recommended to keep it away from sunlight and at a temperature between 10-30 °C.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

Shelf life

2 years.

Packing

Presentation:

25 PCD224-C + 25 BT224 per box.

Packing information: Product description, standards and regulations data, storage conditions, manufacturer information and data on box's label.

Labeling

On product's packing: Product code and description, presentation, standards and regulations data, batch number, batch number of the Biological Indicator, manufacture and expiration dates, bar code and Data Matrix.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Device Industries.

Other relevant information

Read product's Instructions for Use thoroughly before use.

Precautions

Do not store the product near sterilizing agents. Do not use PCDs for monitoring ethylene oxide, dry heat, formaldehyde or any sterilization process other than steam.

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