Infection Control Division | **Chemical Indicators**

KPCD225-C Steam Test Pack PCD Kit

Process Challenge Device for Steam sterilization processes

Usage

For routine load release, especially for loads containing implants.
For routine sterilizer monitoring. For periodic validation of the sterilizer (after installation, repair, relocation).
Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity displacement, steam autoclaves between 132-135 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN, ISO 13485:2016.
ISO 11140-1:2014, ISO 11138-1:2017, ISO 11138-3:2017, ANSI/AAMI ST79.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Classification

Class 1, according to risk (ANMAT).

Characteristics

The kit consists of:

- PCD225-C: Porous cards system holding:
 - A Self-contained Biological Indicator Bionova® BT225.
 - A Type 5 chemical indicator Integron® IT26-C (moving-front integrator).
 - A self-adhesive Record Card, PCDBI-C-RC, where data on sterilization cycle may be written.
- Self-contained Biological Indicator Bionova® BT225 (to be used as a positive control)

Cardboard box (10.6 cm x 6.8 x 1.7 cm) with a Type 1 Process Indicator.

Process Indicator:

Initial color: **blue**.

Final color: **dark grey or black**.

IT26-C Integrator Indicator (within the challenge device):

The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met.

Integration condition is calibrated against the kill time of a 10^6 *Geobacillus stearothermophilus* ATCC®* 7953 spores' population, calculated in BIER (Biological Indicator Evaluator Resistometer).
Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

*ATCC® is a registered trademark of American Type Culture Collection.

Infection Control Division | **Chemical Indicators**

KPCD225-C Steam Test Pack PCD Kit

Process Challenge Device for Steam sterilization processes

BT225 Instant Self-contained Biological Indicator (within the challenge device): $\geq 10^6$ *Geobacillus stearothermophilus* ATCC® 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

Final results: instant fluorescence readout (7 second readout) after incubation at 60 °C (sensitivity $\geq 97\%$). A visual confirmation of the result is optional (color change of the medium caused by pH change) and can be made after incubation at 60 °C for 48 hours and/or for 7 days.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *Geobacillus stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep in a dark place.

Transport conditions

It is recommended to store it away from sunlight and at a temperature between 10 - 30 °C.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

Acceptable excursions during transport are as follows:

T = 2-10 °C and T = 30-37 °C for no longer than seven days. No restrictions for relative humidity.

Shelf life

18 months.

Packing

Presentations:

25 PCD225-C + 25 BT225 Biological Indicators per box.

25 PCD225-C + 5 BT225 Biological Indicators per box.

Packing information: product code and description, regulation, conditions for use, instructions for use, storage conditions, manufacturer information and data on box's label.

Infection Control Division | **Chemical Indicators**

KPCD225-C Steam Test Pack PCD Kit

Process Challenge Device for Steam sterilization processes

Labeling

On product's packing: product code and description, presentation, regulation, batch number, batch number of the Self-contained Biological Indicator, manufacture and expiration date, bar code and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other relevant information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

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