Infection Control Division | Accessories

Bionova® IC10/20 Dual Incubator

For biological indicators and hygiene monitoring systems

Usage

For incubation of:

Bionova® biological indicators: BT10 (EO): 48 h at 37 °C, BT20 (steam): 24 h at 60 °C, BT30 (dry heat): 48 h at 37 °C and BT91 (VH₂O₂): 24 h at 60 °C.

Bionova® self-contained spore ampoules: BT21 (steam): 48 h at 60 °C, BT22 (steam): 48 h at 60 °C, BT23 (steam): 48 h at 60 °C and BT24 (steam): 48 h at 37 °C

Chemdye® hygiene monitoring systems: PRO1 MICRO: 7 minutes at 60 °C and PRO1 ENDO: 7 minutes at 60 °C.

Bionova® Culture medium MC1020-2: intended for use with spore strips and spore coupons biological indicators for VH₂O₂ (24 h at 55–62 °C), steam (24 h at 55–62 °C or 48 h at 60 °C), formaldehyde (48 h at 55–62 °C), EO (48 h at 37 °C) and dry heat (48 h at 37 °C).

Bionova® Culture medium MC1030-2: intended for use with spore strip biological indicators for EO and dry heat (48 h at 37 °C).

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/ EN ISO 13485:2016.

Electrical Safety: IEC 61010-1; IEC 61010-2-010.

Low Voltage Directive 2014/35/EU.

Electromagnetic Compatibility: EN 61326-1; Directive 2014/30/EU.

RoHS Directive 2011/65/EU.

WEEE Directive 2012/19/EU.

Characteristics

Dimensions: 12.8 cm high, 12.4 cm diameter.

10 metal positions (1.3 cm diameter, 3.0 cm deep) to incubate BT23, MC1020-2, MC1030-2, PRO1 ENDO and PRO1 MICRO.

26 metal positions (0.9 cm diameter, 3.0 cm deep) to incubate BT10, BT20, BT30, BT91, BT21, BT22 and BT24.

Special hole for external thermometer.

1 ampoule crusher.

Voltage range: 100 - 240 V AC.

Power: 20 W.

Frequency: 50 - 60 Hz.

Operation temperature: 10 - 30 °C.

Storage conditions

T = 5-40 °C, RH = 30-80 %, keep away from direct sunlight.

Transport conditions

It is recommended to respect the storage conditions during transport.

Products should be transported in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent any risk for human health.

Expected lifetime

5 years.

Packing

1 unit per box.

Packing information: product code and description, regulatory data and manufacturer information.

Infection Control Division | Accessories

Bionova® IC10/20 Dual Incubator

For biological indicators and hygiene monitoring systems

Labelling

On product: product code and description, operating conditions, manufacture date, batch number, serial number, Data Matrix and manufacturer information.

On product's box: product code and description, storage conditions, manufacture date, batch number, serial number and Data Matrix.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Devices industries.

Other important information

Bionova® IC10/20 Incubator has a switch button located at the base of the incubator to select the temperature at $37 \pm 2 \text{ }^\circ\text{C}$ or $60 \pm 2 \text{ }^\circ\text{C}$.

Bionova® IC10/20 incubator allows to calibrate the incubation temperature with the Bionova® TB-IC1020 thermometer. Follow the detailed instructions in the User Manual to calibrate the incubator.

Precautions

Do not pour any liquid inside.

Do not immerse in any liquid.

Only for indoor use.

Disconnect the power cord before cleaning.

Do not use abrasive or corrosive cleaners or disinfectants.

Make sure that the incubator is connected to a properly rated power cord.

Not compatible with Bionova® Rapid, Super rapid, Ultra rapid, Photon nor Hyper, fluorescence readout biological indicators. Use only for extended incubation for an optional visual color change confirmation.

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