Infection Control Division | **Accessories**

TB-IC1020 probe

For temperature sensing

Usage

Accessory for Bionova® Incubators IC10/20, IC10/20FR, IC10/20FRLCD, MiniPro y MiniBio.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

Characteristics

Dimensions: 70 mm (probe).

Range: 0 to + 100 °C / 32 to +212 °F.

Resolution: 0.1 °C/°F.

Sampling Rate: 1 measurement/second.

Battery: 1 pc x 1.5 volt type LR44 or equivalent.

Thermometer size: 125 mm (L).

Environmental operating conditions: 10-30 °C (50-86 °F), RH 30-80 %.

Environmental conditions during manufacture

T = 15-30 °C, RH = 35-85 %.

Storage conditions

T = 10-30 °C, RH = 30-80 %.

Transport conditions

Storage conditions should be strictly followed. The product should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

Shelf life

5 years.

Packing

1 unit per box.

Packing information: product code and description, specifications, measure procedure, considerations, safety information, manufacturer information.

Labeling

On product's packing: product code, calibration certificate, batch number, serial number and manufacturer information.

Possible target markets

Healthcare, Food, Pharmaceutical and medical devices industries.

Infection Control Division | **Accessories**

TB-IC1020 probe

For temperature sensing

Other important information

Use as an accessory for incubators IC10/20, IC10/20FR, IC10/20FRLCD, MiniPro and MiniBio to monitor temperature at 37 ± 2 °C or 60 ± 2 °C.

The metal probe MUST BE completely inserted into the incubator.

Precautions

For indoor use only.

Do not use the thermometer if there are any signs of damage.

Do not bend the probe.

Do not immerse into any liquid

Do not use the thermometer outside its environmental operating conditions.

Do not use the thermometer for continuous temperature monitoring purposes.

Do not store at direct sunlight, high temperature or high humidity.

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