Infection Control Division | Hygiene Indicators

ATP-s1 Surface ATP Test

Hygiene monitoring system based on the detection of residual ATP on surfaces

Usage

Designed for the detection of ATP on hard-to-reach surfaces that have been subjected to cleaning processes.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

Characteristics

- The system has a high absorption swab that allows the collection of samples from different surfaces, and a reactive solution.
- The swab is pre-moistened and is capable of releasing the ATP present inside living cells.
- Product diameter: 12.4 mm \pm 0,5 mm; product length: 205.0 mm \pm 3 mm; swab length: 112.0 mm \pm 3.0 mm.
- Quantitative analysis: using 3M™ Clean-Trace™ Luminometer LM1 and 3M™ Clean-Trace™ Luminometer LX25. Read the instructions for use of the device beforehand. Upon activation of the test, a chemical reaction with the ATP collected on the swab occurs and light is produced. The amount of light produced is proportional to the degree of contamination. Light measurement requires the use of a compatible Luminometer, for example 3M™ Clean-Trace™ Luminometer LM1, results are displayed in Relative Light Units (RLU). The higher the RLU number, the greater the contamination on the analyzed surface.
- This product can be used by personnel with minimal training and is easily integrated into your company's HACCP programs to monitor cleaning effectiveness.

Environmental conditions during manufacture

T = 15-25 °C, RH = 30-80 %.

Storage conditions

T = 2-8 °C, RH = 30-80 %, store in a dark place.

Transport conditions

Storage conditions should be strictly followed.

Transport in closed and strengthened boxes in order to avoid damages. The transport of this product does not represent a risk for health.

Shelf life

12 months when stored refrigerated (2-8 °C) and 1 month when stored at room temperature (\leq 25 °C).

Packing

Presentation: 100 units per box (in 10 bags containing 10 units each).
Packing information: product description, storage conditions, manufacturer information and data on label.

Infection Control Division | Hygiene Indicators

ATP-s1 Surface ATP Test

Hygiene monitoring system based on the detection of residual ATP on surfaces

Labelling

On the product: product description, batch, manufacture date, expiration date and datamatrix code.

On the bag: product code and description, storage conditions, manufacturer information, presentation, batch number, manufacture date, expiration date, bar code and datamatrix code.

On product's box: Product code, batch number, manufacture and expiration date, barcode and datamatrix code.

Possible target markets

Healthcare, food, dental, veterinary, pharmaceutical and medical industries.

Other important information

Read product's instructions for use thoroughly before use. The device is compatible with 3M™ Clean-Trace™ Luminometer.

Precautions

Do not open until its use.

Do not reuse.

Do not touch the swab with your hands. Wear gloves, surgical mask and other hygiene measures.

Do not freeze.

High concentrations of cleaning products or disinfectants can affect the development of the reaction.

Keep protected from light until use.

Let it cool to room temperature before use.

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