Infection Control Division | **Cleaning Indicators**

## *splat* CDWA4 Cleaning Indicator

For cleaning performance tests.

### Usage

Designed for routine control of cleaning processes and for monitoring automatic washing machines effectiveness as well as monitoring ultrasonic washing machines which operate at frequencies equal to or greater than 35 kHz, and temperatures between 30 °C and 70 °C.

### Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.  
ISO/TS 15883-5:2005.

### Characteristics

Challenge level: very high.

Allows detection of all factors affecting the outcome of cleaning process.

Initial Color: **fuchsia**.

Size: 70 mm x 30 mm (approximate).

Synthetic temperature-stable support.

Result: elimination of fuchsia organic substance.

To interpret results visually, use the results reference guide on product's instructions for use.

Optionally, the readout of the cleaning indicator can be done using Trazanto Automatic Quality Control and Traceability System (see instructions for use in device's manual).

**NOTE:** If unfavorable results are obtained, check quality/quantity of detergent, cleaning time and temperature, positioning of the load, cleaning indicator location and general condition of the washing machine.

### Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

### Storage conditions

T = 10-30 °C, RH = 30-80 %, keep in a dark place in its original box.

### Transport conditions

It is recommended to keep it away from sunlight and at a temperature between 10 - 30 °C.

Transport in closed and strengthened boxes in order to avoid damages.

Transportation of this product does not represent a risk for health.

Infection Control Division | **Cleaning Indicators**

## *splat* CDWA4 Cleaning Indicator

For cleaning performance tests.

### **Shelf life**

5 years.

### **Packing**

200 units per aluminum foil bag.

Packing information: product description, intended use, instructions for use, standards reference, storage conditions, manufacturer information.

### **Labelling**

On product's label: product's code and description, product's presentation, batch number, expiration date, bar code and datamatrix code.

### **Possible target markets**

Healthcare, Food, Pharmaceutical and Medical industries.

### **Other important information**

Read product's instructions for use thoroughly before use.

### **Precautions**

Store the product within its' original package, away from disinfectants, corrosive or oxidant agents, gases and humidity.

Do not re-use the product.

Do not touch indicator ink with your hands.

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