CDWU

Ultrasonic Cleaning Indicator

For ultrasonic cavitation performance tests.

Usage

Designed for ultrasonic cleaning processes' routine control and monitoring ultrasonic cavitation effectiveness of ultrasonic washing machines operating at frequencies equal to or greater than 35 kHz, at temperatures ranging from 18 °C to 70 °C.

Authorization

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

Characteristics

2 ml clear glass vial: 31.0 mm high x 11.6 mm external diameter, 8.0 mm neck diameter. Wall thickness: 0.5 mm.

Brown/Red rubber septum: 8.0 x 1.0 mm.

Aluminum top cap: 5.5 mm high x 11.7 mm diameter. Wall thickness: 0.3 mm.

Glass beads: 3.0 mm.

Blue clear reactive solution: 1.5 ml.

Allows detection of deficiencies in the cavitation performance of the ultrasonic washing machine.

A self-adhesive Record Card accompanies each vial where data about the ultrasonic washing control process can be recorded.

Challenge level: high

Initial Color: blue

Result: yellow transparent solution

NOTE: If unfavorable results are obtained, check working frequency of the ultrasonic cleaner, time and temperature of the ultrasonic washing, presence of bubbles within the wash solution, location of indicator and general condition of ultrasonic washing machine.

Environmental conditions during manufacture

T = 10-30 °C, RH = 30-80%.

Storage conditions

T = 10-30 °C, RH = 30-80%, store in dark and dry place inside original packing.

Transportation conditions

Storage conditions should be strictly followed.

Transport in closed and strengthened boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

18 months.

Packing

6 or 30 units per box.

Packing information: product code and description, manufacturer information and data on label.

Labelling

On product's box: product code and description, batch number, manufacture and expiration date, barcode and datamatrix code.

Possible target markets

Healthcare and Industry.

Other important information

Read product's directions for use thoroughly before use.

Description of identified and/or predictable risks which could represent a risk for the product

Do not re-use the CDWU indicator.

Do not open the indicator before its use, since this will cause its deterioration.

Do not use the CDWU indicator for monitoring normal washing using non-ultrasonic cleaners.

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