Infection Control Division | **Cleaning Indicators**

## KPRO2-E69 Hygiene Monitoring System

Kit system for detecting residual protein on surfaces

### Usage

Designed to detect proteins on surfaces of surgical equipment after the cleaning process and in hard to reach areas.


### Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 15883-1:2006 / HTM 01-01:2016 / HTM 01-05:2013

### Characteristics

- The system has a high absorption 6.9 cm long swab, SW69, as well as a plastic tube with reactive solution. The SW69 swab allows the collection of samples from different surfaces with the same efficacy following the cleaning process.
- Volume of reactive solution in the plastic tube: 0.45 ml.
- An instant and visual reading of a color change indicates the presence of detectable levels of protein. When comparing the color obtained in the solution with the reference color card, a visual estimation of surface cleaning can be made.
- KPRO2-E69 can detect as little as 1 µg of protein.

### Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

### Storage conditions

T = 2-30 °C, RH = 30-80 %, keep in a dark place.

### Transport conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent any risk for human health. Product can be kept at temperature below 25 °C for up to 25 days during transport without modifying its shelf-life (provided that the subsequent storage conditions are respected).

### Shelf life

24 months when stored refrigerated (2 °C - 8 °C) and 6 months when stored at room temperature (up to 30 °C).

Infection Control Division | **Cleaning Indicators**

## KPRO2-E69 Hygiene Monitoring System

Kit system for detecting residual protein on surfaces

### Packing

18 plastic tubes with reactive solution + 18 SW69 swabs + 2 moisturizer tubes for swabs + 18 registration cards + 1 reference color card.

Packing information: product code and description, storage conditions, manufacturer information and data on label.

### Labelling

On swabs bag: product code and description, batch number and presentation.

On product's box: Product code and description, batch number, manufacture and expiration date, presentation, barcode and datamatrix code.

### Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

### Other important information

Read product's instructions for use thoroughly before use.

### Precautions

Do not open until its use.

Do not reuse.

Do not touch the swab with your hands. Wear gloves, surgical mask and other hygiene measures.

Do not freeze.

Although KPRO2 has high sensitivity, it is not a microbiological monitoring device, and surfaces may not be considered to be sterile based on the results obtained.

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