Infection Control Division | Hygiene Indicators

PRO1 VT Validation Kit

System for efficacy validation of the PRO1 MICRO Hygiene Monitoring System

Usage

Designed for efficacy validation of Chemdye® PRO1 MICRO Hygiene Monitoring System when used with Bionova® MiniPro, IC10/20FR or IC10/20FRLCD Auto-Reader Incubators. The PRO1 MICRO allows to quantify with high sensitivity proteins and reducing agents on different surfaces and difficult access areas after the cleaning process.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.
ISO 15883-1: 2006 / HTM 01-01: 2016 / HTM 01-05: 2013.

Characteristics

The system has a Negative and a Positive Control.

The Negative Control consists of a vial with a protein-free solution (transparent cap).

The Positive Control consists of a vial with a defined amount of standard protein (light blue cap).

Vial's diameter: 12.0 mm; Vial's length: 47.0 mm.

Quantitative analysis: Using Bionova® IC10/20FR, IC10/20FRLCD or MiniPro Auto-reader incubators. Read the instructions for use of the device previously.

Incubation time: using Bionova® IC10/20FR or Bionova® IC10/20FRLCD: 7 minutes; using Bionova® MiniPro: 4 minutes.

Incubation range: $(60 \pm 2) ^\circ\text{C}$.

Environmental conditions during manufacture

15 - 30 °C, RH 30 - 80 %.

Storage conditions

T = 2-8 °C, RH 30 - 80 %, store in a dark place.

Transport conditions

It is recommended to store it away from sunlight and at a temperature between 2 - 8 °C.

Transport in closed and strengthened boxes in order to avoid damages. The transport of this product does not represent a risk for health.

Acceptable excursions during transport are as follows:

Room temperature (25 - 30 °C) for no longer than fifteen days. No restrictions for relative humidity.

Shelf life

Shelf life of 2 years when stored refrigerated (between 2 - 8 °C).

Infection Control Division | Hygiene Indicators

PRO1 VT Validation Kit

System for efficacy validation of the PRO1 MICRO Hygiene Monitoring System

Packing

Box with 10 Negative Controls (transparent cap) and 10 Positive Controls (light blue cap).

Packing information: product code and description, process for intended use, presentation, storage conditions, manufacturer information and data on label.

Labelling

On product's pack: product code and description, batch number, manufacture date and expiration date, barcode and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information

Read product's instructions for use thoroughly before using.

Precautions

Do not reuse.

Do not open until its use.

Do not touch the swab with your hands. Wear gloves, surgical mask and other hygiene measures.

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