Infection Control Division | Hygiene Indicators

PRO1 ENDO Hygiene Monitoring System

High sensitivity pen system for the detection and quantification of proteins in cannulated instruments

Usage

Chemdye® PRO1 ENDO Hygiene Monitoring System is designed to verify the cleanliness of cannulated instruments by detecting protein residues that may remain after the cleaning process. The system is used in conjunction with Chemdye® SWE high absorption swabs, specially designed for cannulated instruments.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

Characteristics

The system consists of a reactive solution and is used in conjunction with a high absorption swab, which allows the collection of samples from the internal channel of cannulated instruments.

- Limit of quantification: 1.0 µg of BSA.
- Quantitative analysis: using the Bionova® IC10/20FR, Bionova® IC10/20FRLCD or Bionova® MiniPro auto-reader incubators. Read the instructions for use of the device previously.

Incubation time: using the Bionova® MiniPro auto-reader incubators: 4 minutes; using the Bionova® IC10/20FR or the Bionova® IC10/20FRLCD auto-reader incubators: 7 minutes.

Incubation range: 60 ± 2 °C.

- Qualitative analysis: incubate using the Bionova® IC10/20 incubator (temperature: 60 ± 2 °C; 7 minutes) and interpret visually by using the reference guide available inside the product's box.

Storage conditions

Chemdye® PRO1 ENDO System should be stored at temperatures between 2-8 °C, protected from light and relative humidity (RH) of 30-80 %. Store vertically as indicated on the packaging.

Transport conditions

It is recommend to respect storage conditions during transport.

Transport in closed and strengthened boxes in order to avoid damages. The transport of this product does not represent a risk for health.

Shelf life

24 months when refrigerated (2-8 °C). It can also be stored at room temperature (25-30 °C); however, its shelf life is reduced to 6 months in this case.

Packing

Box with 20 units (1 aluminum bag with 20 pens + 4 vials of moisturizer solution).

Packing information: product code and description, presentation, storage conditions and manufacturer information.

Labelling

On product: product code and description, batch number, manufacture and expiration date and datamatrix code.

On product's box: product code, batch number, manufacture and expiration date, bar code and datamatrix code.

Infection Control Division | Hygiene Indicators

PRO1 ENDO Hygiene Monitoring System

High sensitivity pen system for the detection and quantification of proteins in cannulated instruments

Possible target markets

Healthcare, Pharmaceutical and Medical industries.

Other important information

Read product's instructions for use thoroughly before use.

Precautions

Handle with gloves at all times.

Do not open until use.

Do not reuse. Do not re-incubate.

Do not swab surfaces that have not been subjected to a cleaning process or visibly contaminated areas.

Avoid touching the absorbent part of the swab or the inside of the device with your hands. Use gloves, a mask, and other hygiene measures.

Do not freeze.

Do not use Chemdye® SWE swabs on any device other than Chemdye® PRO1 ENDO.

High concentrations of alkaline detergents may cause false-negative results in some cases. Oxidizing disinfectants (such as hydrogen peroxide or peracetic acid) may cause a false-positive color change reaction.

Although Chemdye® PRO1 ENDO is highly sensitive, its intended use is for monitoring channels of cannulated instruments after the cleaning process, therefore the test result is not evidence that the device is sterile.

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