# CD202

## Multivariable Chemical Indicator

*For monitoring Steam sterilization processes*

### Usage

Control of Steam sterilization processes.

**Conditions:** 121 °C, 15 minutes // 134 °C, 3.5 minutes.

### Applicable regulation

Designed under Quality Management Systems standards ISO 13485:2016/NS-EN ISO 13485:2016. ISO 11140-1:2014

### Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-4.

### Classification

Class 1, according to risk (ANMAT).

### Characteristics

Type 4 Multivariable Chemical Indicator.


Punched and one-side-laminated double card strips printed with indicator ink.

Size: 200 x 15 mm

Initial color: **blue**

Final color: **black**

100 % Toxic Heavy Metals free

### Environmental conditions during manufacture

T = 15-30 °C, RH 30-80 %.

### Storage conditions

T = 10-30 °C, RH 30-80 %, keep out of direct light.

### Transport conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent a risk for health.

### Shelf life

5 years

### Packing

250 double strips per aluminum foil bag.

Packing information: product description, instructions for use, regulation, storage conditions, manufacturer information and data on label.

### Labelling

On product's packing: product code and description, process for intended use, classification according to regulation, color change, presentation, batch number, manufacture and expiration date, bar code and datamatrix code.

**Possible target markets**

Healthcare, Food, Pharmaceutical and Medical industries.

**Other important information**

Read product's instructions for use thoroughly before use.

**Precautions**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat or any sterilization processes other than Steam.

Do not reuse the sterilizer until the indicator has changed into the color indicated in the product.

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