Infection Control Division | Chemical Indicators

CD42 Process Indicator

For monitoring Plasma or Vaporized Hydrogen Peroxide sterilization processes

Usage

United States

Terragene® Chemdye® CD42 is a Chemical Process Indicator intended for monitoring the efficacy of Vaporized Hydrogen Peroxide sterilization processes. The chemical indicator changes from red to yellow to indicate that the conditions of the cycle have been met. Intended use cycles: STERRAD 100S 54 minutes, STERRAD 100NX Standard and Express cycles and V-Pro Max and Sterizone VP4 cycles.

Outside the United States

Control of plasma or vaporized H₂O₂ sterilization processes.

Conditions: 6 minutes, 50 °C, 2,3 mg/l H₂O₂.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN, ISO 13485:2016.

ISO 11140-1:2014.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Classification

Class 1, according to risk (ANMAT).

FDA 510(k)

K191021.

Characteristics

Type 1 Chemical Indicator (according to ISO 11140-1).

Varnished strips printed with indicator ink.

Size: 105 x 18 mm (approximate).

Initial color: **Red.**

Final color: **Yellow.**

100 % Toxic Heavy Metal free

To interpret results visually, use the reference guide found in the Indications for Use.

Optionally, the readout of the chemical indicator can be done using Trazanto Automatic Quality Control and Traceability System (see instructions for use in device's manual).

Infection Control Division | Chemical Indicators

CD42 Process Indicator

For monitoring Plasma or Vaporized Hydrogen Peroxide sterilization processes

Environmental conditions during manufacture

T= 15-30 °C, RH= 30-80 %.

Storage conditions

T= 10-30 °C, RH= 30-80 %, keep in a dark place in its original box.

Transport conditions

It is recommended to store it away from sunlight and at a temperature between 10-30 °C.

Transport in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent a risk for health.

Acceptable excursions during transport are as follows:

T = 2-10 °C and T = 30-37 °C for no longer than seven days. No restrictions for relative humidity.

Shelf life

5 years.

Packing

250 or 500 strips per aluminum foil bag.

Packing information: product description, instructions for use, storage conditions and manufacturer information.

Labelling

On product's packing: product code and description, process for intended use, classification according to regulation, color change, product presentation, batch number, expiration date, barcode and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to Steam, Dry Heat or any sterilization processes other than H₂O₂.

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