Infection Control Division | Chemical Indicators

IT26-C Moving front chemical integrator

For Steam sterilization processes

Usage

United States

Chemically react over time with the critical parameters of steam sterilization cycle within a specified tolerance. Intended use cycles: Gravity Displacement: 121 °C for 30 minutes, 132 °C for 15 minutes, 132 °C for 25 minutes, 135 °C for 10 minutes; Dynamic Air Removal (Vacuum Assist): 132 °C for 4 minutes, 135 °C for 3 minutes.

SV121 °C/16.5 min., SV132 °C/2.0 min., SV135 °C/1.2 min.

Outside the United States

Control of Steam sterilization processes between 118 °C and 138 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN, ISO 13485:2016.

ISO 11140-1:2014.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Classification

Class 1, according to risk (ANMAT).

FDA 510(k)

K191021.

Characteristics

Type 5 Integrating Indicator (according to ISO 11140-1).

Size: 52 mm x 21 mm.

Chemical Integrator consists of a paper wick and a steam and temperature sensitive chemical pellet contained in a paper/film/foil laminate. The chemical pellet melts and migrates as a dark bar along the paper wick. The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met. The extent of migration depends on steam quality, time, and temperature. The ACCEPT result is reached when a theoretical spore population reaches its kill time, indicating integration condition has been reached. This condition is calibrated with the kill time of 10^6 *Geobacillus stearothermophilus* ATCC 7953* spore population, calculated in BIER (Biological Indicator Evaluator Resistometer).

*ATCC® is a registered trademark of American Type Culture Collection.

Infection Control Division | **Chemical Indicators**

IT26-C Moving front chemical integrator

For Steam sterilization processes

Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

Optionally, the readout of the chemical indicator can be done using Trazanto Automatic Quality Control and Traceability System (see instructions for use in device´s manual).

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep out of direct light.

Transport conditions

It is recommended to store it away from sunlight and at a temperature between 10 - 30 °C.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

Acceptable excursions during transport are as follows:

T = 2-10 °C and T = 30-37 °C for no longer than seven days. No restrictions for relative humidity.

Shelf life

5 years.

Packing

250 units per aluminum foil bag.

Labelling

On product´s packing: product code and description, process for intended use, instructions for use, storage conditions, product presentation, classification according to regulation, manufacturer information, batch number, manufacture and expiration dates, SV at 121 °C, 128 °C and 135 °C, barcode and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information

Read product´s instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents, strong alkaline or acidic products such as cleaning/disinfecting agents.

Do not expose this product to EO, Dry Heat or any sterilization processes other than Steam.

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