IT26SAD

Unique point Integrator Indicator with self-adhesive support

For Steam sterilization processes

Usage

Control of Steam sterilization processes between 121 °C and 135 °C

Applicable regulation

Designed under Quality Management Systems standards ISO 13485:2016/NS-EN ISO 13485:2016. ISO 11140-1:2014

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-4.

Classification

Class 1, according to risk (ANMAT).

Characteristics

Type 5 Integrator Indicator.

Strips printed with indicator ink. Version of the integrator indicator in a self-adhesive support, which allows their easy record.

Size: 72 mm x 23 mm

Initial color: **purple**

Final color: **green**

Integration condition is calibrated with the kill time of a 10^6 *Geobacillus stearothermophilus* spores population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

The product ensures an adequate control of the efficacy of sterilization processes (temperature, time, quality of steam).

100 % Toxic Heavy Metals free

Environmental conditions during manufacture

T = 15-30 °C, RH 30-80 %.

Storage conditions

T = 10-30°C, RH 30-80 %, keep out of direct light.

Transport conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent a risk for health.

Shelf life

5 years

Packing

100 units per aluminium foil bag.

Packing information: product description, regulation, storage conditions, manufacturer information and data on label.

Labelling

On product's packing: product code and description, process for intended use, classification according to regulation, color change, presentation, batch number, manufacture and expiration date, bar code and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat or any sterilization processes other than Steam.

Do not reuse the sterilizer until the indicator has changed into the color indicated in the product.

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