Infection Control Division | **Chemical Indicators**

KBD8948X Bowie-Dick Test Card Kit

Includes Bowie-Dick Test Card Chemical Indicator and its Holder

Usage

Bowie-Dick test monitors the effectiveness of air removal in vacuum-assisted steam sterilizers. Designed for cycles of 3.5 minutes at 134 °C and 4 minutes at 132 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN, ISO 13485:2016, ISO 11140-1:2014 and ISO 11140-4:2007.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Classification

Class 1, according to risk (ANMAT).

Characteristics

BD8948X

Bowie-Dick Test Card Chemical Indicator.

Type 2 Chemical Indicator.

Initial Color: **purple.**

Final Color: **green.**

Size: 89 mm x 48 mm x 0.33 mm

Material: one-side-laminated cardboard.

Cycle information can be recorded on the back of the card.

NOTE: Air removal problems are shown as heterogeneous color change, i.e. the presence of purple color.

100 % Toxic Heavy Metal free

BD8948H

Stainless steel re-usable Holder for Bowie-Dick Test Card of great resistance and durability coated by a protective layer.

Keeps BD8948X Test Card in place for proper assessment of sterilization cycle.

Height: 42.7 mm

Base: 133.2 mm x 72.4 mm

Upper part where the BD8948X Test Card fits: 93 mm x 53.5 mm.


Environmental conditions during manufacture

T= 15-30 °C, RH= 30-80 %.

Infection Control Division | **Chemical Indicators**

KBD8948X Bowie-Dick Test Card Kit

Includes Bowie-Dick Test Card Chemical Indicator and its Holder

Storage conditions

T= 10-30 °C, RH= 30-80 %, keep out of direct light.

Transport conditions

It is recommended to store it away from sunlight and at a temperature between 10-30 °C.

Transport in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent a risk for health.

Acceptable excursions during transport are as follows:

T = 2-10 °C and T = 30-37 °C for no longer than seven days.

No restrictions for relative humidity.

Shelf life

5 years.

Packing

1 Holder + 120 Bowie-Dick Test Cards per box.

Packing information: product description, process for intended use, manufacturer information, storage conditions and data on the label.

Labelling

On product's packing: Product code, classification according to regulation, color change, presentation, batch number, manufacture date, expiration date, barcode and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat or any sterilization process other than Steam.

Do not use a Bowie-Dick Test Card in gravity displacement steam sterilization cycles (where there is no pre-vacuum phase) or any kind of sterilization processes other than vacuum-assisted Steam.

Do not use the sterilizer until the indicator has reached its homogeneous endpoint color.

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