
# KH2X15-9.0/M

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
# KH2X15-9.0/M

## Helix-PCD + Chemical Indicator test strips Kit.

*For Steam sterilization processes control.*

### Usage

Process Challenge Device for Steam sterilization processes control.  
For safe release of hollow loads.  
Conditions: 9 min at 134 °C

### Applicable regulation

EN 867-5, EN 285, ISO 11140-1:2014

### Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012.

### Characteristics

Helix-PCD

Metal capsule of great resistance and durability, to hold a Chemical Indicator designed for this system. The capsule is connected to a high resistance plastic tube of 1.5 m length (see catalogue for other lengths) and 2 mm internal diameter. Special silicone connector coated by external connector.

PCD-A-9.0

Type 2 Chemical Indicator.

Self-adhesive indicator ink printed strip. Size: 77 x 6 mm.

10 strips per sheet. Size: 80 x 100 mm.

Initial color: purple

Final color: green

NOTE: If the 4 bars change to the final color, the sterilization cycle was successful, the load can be securely released. If one or more bars did not turn to the final color, this shows a failure in sterilizing agent penetration, the load cannot be considered sterile, that is, it cannot be released.

100 % Toxic Heavy Metals free

### Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

### Storage conditions

T = 10-30 °C, RH 30-80 %, store in dark place.

### Transportation conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

### Shelf-life

Metal capsule: 500 cycles. Plastic helix: 250 cycles (replaceable at very low cost). PCD-A-9.0: 5 years

**Packing**

250 PCD-A-9.0 strips in sheets of 10 units + 1 Helix-PCD + 1 cotton bag.

Packing information: product code and description, process for intended use, presentation, manufacturer information and data on label.

**Labelling**

On product's box: product code and batch number.

On helix-PCD bag: product code and batch number.

On PCD-A-9.0 box: batch number and expiration date.

**Possible target markets**

Healthcare, Food, Pharmaceutical and Medical industries.

**Other important information (optional)**

Read product's directions for use thoroughly before use.

**Note:** When necessary, limits of each product's parameters should be specified.

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Formaldehyde or any sterilization processes other than Steam.

Do not reuse the sterilizer until the color change of the chemical indicator is the one indicated in the product.