# PCD-A-7.0Y

## Chemical Indicator for use in Helix-PCD system

*For Steam sterilization processes control.*

### Usage

Designed to be used in the Helix-PCD system for secure load release.

**Conditions:** 7.0 minutes, 134 °C // 20 minutes, 121 °C.

### Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11140-1:2014.

### Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM-1614-4.

### Classification

Class 1, according to risk (ANMAT).

### Characteristics

Type 2 Chemical Indicator.

Self-adhesive varnished strip printed with indicator ink.

Size: 76 x 6 mm.

10 strips per sheet. Size: 80 x 98.5 mm.

Initial color: **yellow**

Final color: **dark brown or black**

100 % Toxic Heavy Metal free

NOTE: If the 4 bars change to the final color, the sterilization cycle was successful, the load can be securely released. If one or more bars did not turn to the final color, this shows a failure in sterilizing agent penetration, the load cannot be considered sterile, that is, it cannot be released.

### Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

### Storage conditions

T = 10-30 °C, RH = 30-80 %, keep out of direct light.

### Transport conditions

Storage conditions should be strictly followed. Transport in closed and strengthened boxes in order to avoid damages. The transport of this product does not represent a risk for health.

### Shelf life

5 years

### Packing

250 strips in sheets of 10 units.

Packing information: product description, classification according to regulation, storage conditions, manufacturer information and data on label.

**Labelling**

On product packing: product code, process for intended use, presentation, color change, batch number, manufacture and expiration date, barcode and datamatrix code.

**Possible target markets**

Healthcare, Food, Pharmaceutical and Medical industries.

**Other important information**

Read product's instructions for use thoroughly before use.

**Precautions**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Formaldehyde or any sterilization process other than Steam.

Do not reuse the sterilizer until the color change of the chemical indicator is the one indicated in the product.

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