PCD-A-3.5BD

Bowie-Dick simulation test for use in Helix-PCD system

For Steam sterilization processes control.

Usage

Designed to be used in the Helix-PCD system for assessing air removal and steam penetration of vacuum-assisted steam sterilizers.

Conditions: Bowie-Dick cycle: 3.5 minutes, 134 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11140-1:2014.

Product Classification

Class 1, according to risk.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM-1614-4.

Characteristics

Type 2 Chemical Indicator.

Self-adhesive strip printed with indicator ink. Size: 76 x 6 mm.

10 strips per sheet. Size: 80 x 98.5 mm.

Initial color: **blue**

Final color: **dark grey/black**

100 % Toxic Heavy Metal free

Note: If 4 bars change to the final color, the result is successful, the sterilizer works properly. If one or more bars have not turned to the final color, this shows a failure in sterilizing agent penetration, therefore an equipment malfunction.

Environmental conditions during manufacture

T = 15-30 °C, RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep out of direct light.

Transport conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent a risk for health.

Shelf life

5 years

Packing

250 strips in sheets of 10 units.

Packing information: product description, classification according to regulation, storage conditions, manufacturer information and data on label.

Labelling

On product packing: product code, process for intended use, presentation, color change, batch number, manufacture and expiration date, barcode and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Formaldehyde or any sterilization process other than Steam.

Do not reuse the sterilizer until the color change of the chemical indicator is the one indicated in the product.

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