PCD-A-E1

Test strip with chemical indicator

For use in Helix-PCD system.

Usage

Chemical Indicators designed to be used in the Helix-PCD system for control of Ethylene Oxide (EO) sterilization processes. **Condition:** 240 min., 54 °C, EO 600 mg/l, 60 % RH.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.
ISO 11140-1:2014.

Product Classification

Class 1, according to risk.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM-1614-4.

Characteristics

Type 2 Chemical Indicator.

Self-adhesive varnished strip printed with indicator ink.

Size: 76 x 6 mm.

10 strips per sheet. Size: 80 x 98.5 mm.

Initial color: **light blue**

Final color: **green**

100 % Toxic Heavy Metals free

NOTE: If the 4 bars change to the final color, the result is successful, the sterilizer is working properly. If any bar color is different from the final color, it shows an equipment malfunction.

Environmental conditions during manufacture

T = 15-30 °C, RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep out of direct light.

Transport conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent a risk for health.

Shelf life

5 years

Packing

250 strips in sheets of 10 units.

Packing information: product description, classification according to regulation, storage conditions, manufacturer information and data on label.

Labelling

On product packing: product code, process for intended use, presentation, color change, batch number, manufacture and expiration date, barcode and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to Steam, Dry Heat, Formaldehyde or any sterilization process other than EO.

Do not reuse the sterilizer until the color change of the chemical indicator is the one indicated in the product.

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