Infection Control Division | **Chemical Indicators**

PCD26-C Process Challenge Device

Test pack for Steam sterilization processes.

Usage

Monitoring dynamic air removal and gravity displacement Steam sterilization processes between 121-135 °C.

Applicable regulation

Designed under Quality Management Systems standards ISO 13485:2016/NS-EN, ISO 13485:2016, ISO 11140-1:2014. ANSI/AAMI ST79.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Classification

Class 1, according to risk.

FDA 510(k)

K191021.

Characteristics

Porous card system holding an IT26-C Chemical Integrator (Type 5 according to ISO 11140-1: 2014). The system also contains a self-adhesive Record Card, PCD26-C-RC, where data on sterilization cycle may be registered.

Cardboard box (10.5 cm x 6.5 cm x 1.8 cm) with a self-adhesive Type 1 Process Indicator.

Process Indicator:

Initial color: **blue**.

Final color: **dark grey or black**.

IT26-C Chemical Integrator (within the challenge device):

The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met.

Integrating condition is calibrated against the kill time of 10^6 *Geobacillus stearothermophilus* ATCC 7953®* spore population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C and 135 °C.

100 % Toxic Heavy Metal free.

*ATCC® is a registered trademark of American Type Culture Collection.

Infection Control Division | **Chemical Indicators**

PCD26-C Process Challenge Device

Test pack for Steam sterilization processes.

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep out of direct light.

Transport conditions

It is recommended to respect the storage conditions during transport.

Transport in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent a risk for health.

Shelf life

5 years.

Packing

Packing information: product code and description, regulation, process for intended use, conditions and instructions for use, storage conditions and manufacturer information.

Labeling

On the packing: product code, batch number of the PCD, batch number of Chemical Integrator, manufacture date, expiration date, barcode and data matrix code.

On PCD26-C-RC: product code, classification of the chemical integrator according to regulation and process for intended use.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other relevant information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat or any sterilization processes other than Steam.

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