

# 24-Month shelf life of Self-Contained Biological Indicators (SCBIs)

## ABSTRACT

This technical note outlines the scientific, regulatory, and economic rationale justifying the universal adoption of a 24-month shelf life for self-contained biological indicators (SCBIs) based on *Geobacillus stearothermophilus* and *Bacillus atrophaeus*. It establishes that while bacterial spores possess an intrinsic viability exceeding 5 to 10 years, the primary limiting factor is the stability of the integrated culture medium. Beyond this timeframe, the medium is susceptible to photo-oxidation of pH indicators, Maillard reactions in nutrients, and pH drift due to glass leaching. Additionally, the physical aging of the multi-material system—comprising the plastic body, adhesives, filtration membranes, and glass ampoules—is accounted for, ensuring all interfaces remain within specifications under real-world conditions. Ultimately, extending shelf life to 36 months requires an additional year of real-time stability studies that offers no proportional commercial return due to typical hospital inventory turnover, rendering the 24-month claim the optimal economic balance and the established benchmark for regulatory convergence.

## 1. Purpose

This technical note explains the scientific, regulatory, and economic rationale behind the 24-month shelf life universally adopted for self-contained biological indicators (SCBIs) used in sterilization monitoring. It addresses the main degradation mechanisms of each system component, compares alternative shelf life claims, and summarizes the convergent factors that make 24 months the optimal declaration for this product category.

## 2. Scope

Applicable to all SCBIs based on *Geobacillus stearothermophilus* (steam sterilization, 121°C / 134°C) and *Bacillus atrophaeus* (ethylene oxide sterilization), incorporating a spore carrier strip and a self-contained culture medium ampoule within a single polypropylene or polycarbonate body. The considerations in this note apply to products from Terragene® Bionova®, 3M Attest, Mesa Labs SGM, and equivalent formats from other manufacturers.

## 3. Background

The 24-month shelf life observed almost universally in SCBIs is neither an arbitrary value nor a regulatory standard imposed by norms, but rather the convergent result of several technical, commercial, and regulatory constraints that converge at that time horizon. Understanding each individual factor also explains why some products reach only 18 months, and why extending to 36 months requires disproportionately greater experimental effort.

## 4. Technical factors governing shelf life

### 4.1 Primary limiting factor: culture medium stability

A properly processed bacterial spore of *G. stearothermophilus* or *B. atrophaeus* — sporulated under optimal conditions, washed, lyophilized or adsorbed onto Whatman No. 1 filter paper, and stored dry at room temperature — has an intrinsic shelf life that easily exceeds 5 to 10 years without significant population loss or detectable change in D-value. Abundant literature (Russell, Setlow, and manufacturers' own historical internal studies) documents viable *Bacillus* spores after decades of storage under appropriate conditions. The spore itself is therefore not the bottleneck for SCBI shelf life.

The true limiting factor is the culture medium embedded within the self-contained vial, particularly in sealed glass ampoule or plastic reservoir formats. The critical quality attributes (CQAs) of the medium that degrade most rapidly are:

Colorimetric pH indicator (bromocresol purple, phenol red, or proprietary combinations): undergoes photo-oxidation and slow chemical degradation even when protected from light.

Nutrient stability (soy peptone, yeast extract, glucose): subject to Maillard reactions, amino acid oxidation, and slow caramelization.

Medium pH: drifts due to buffer hydrolysis and ionic migration from the glass (typical alkaline leaching from low-quality Type I or Type III glass).

Growth promotion capability: erodes as critical nutrients for spore germination (L-alanine, inosine, available manganese) degrade or become sequestered.

Real-time stability studies consistently show that between 24 and 36 months, a progressive

increase in positive control outgrowth time begins to appear, along with less distinct or abnormal colorimetric changes and eventually growth failures in a non-negligible fraction of replicates. Twenty-four months represent the point at which the manufacturer can declare shelf life with statistical confidence (95% CI per ICH Q1E) and without risk of nonconformity under storage conditions reasonably deviated from optimum (global transportation, tropical climate distribution, periods on customer shelves without strict environmental control).

#### **4.2 Secondary factor: integrity of the multi-material system**

SCBIs integrate within a single plastic body a spore carrier paper strip and a glass ampoule containing culture medium. The system contains multiple interfaces that age over time:

Plastic body (polypropylene / polycarbonate): undergoes stress relaxation and may subtly deform, altering the ampoule crushing geometry. Label adhesive: undergoes cold flow and may compromise labeling integrity.

Hydrophobic filtration membrane cap (allows steam or EO ingress during sterilization): may accumulate residual moisture or contaminants.

Internal glass ampoule: may develop microfractures due to cyclic thermal stress during transportation.

Each of these interfaces has its own degradation kinetics, and 24 months represents the envelope within which all of them remain demonstrably within specification under real-world use conditions, not only in controlled warehouse storage.

### **5. Economic-experimental factor: marginal cost of shelf life extension**

Establishing shelf life requires real-time studies that necessarily take as long as the declared duration. To claim 24 months, 24 months of data on three lots are required; to claim 36 months, 36 months of data are required. Extending from 24 to 36 months delays product launch by an additional 12 months — or requires maintaining three additional lots in continuous stability studies for another year before updating the dossier and labeling — with significant opportunity cost and without proportional commercial return, given that hospital inventory turnover rarely exceeds 12 to 18 months in clinical practice. Twenty-four

months represent the optimal cost-benefit balance between development effort, regulatory risk, and real-world utility for the end user.

### **6. Regulatory factor: convergence with established practice**

Although neither ISO 11138-1 nor the FDA impose a specific shelf life, the historical market convergence around 24 months — established by 3M Attest since the 1980s and followed by virtually all manufacturers — created a de facto standard that regulatory authorities use as a comparative benchmark. A manufacturer declaring 12 months would be perceived as inferior; one declaring 60 months would face significantly greater regulatory scrutiny and higher market skepticism, requiring disproportionately extensive evidence packages. Twenty-four months are simultaneously the market expectation and the threshold of reasonable evidence for regulatory justification without friction.

### **7. Known exceptions and comparative analysis**

The following table summarizes BI formats with shelf lives diverging from the 24-month standard, confirming that the culture medium — not the spore — is the determining factor.

### **8. Summary**

The 24-month shelf life of an SCBI is the convergent result of four independent factors:

- Technical: Real chemical stability of the culture medium — the primary limiting factor. Degradation of the colorimetric indicator, nutrient components, and growth promotion capability systematically constrains the useful life of the integrated system.
- Robustness: Integrity envelope of the multi-material system. All degradation interfaces (plastic body, adhesive, membrane, glass ampoule) remain demonstrably within specification within the 24-month window under real-world use conditions.
- Economic: Optimal cost-benefit balance between development time and commercial utility. Extending to 36 months requires 12 additional months of real-time studies with negligible commercial return given typical hospital inventory cycles.
- Regulatory/commercial: Convergence with

market practice and regulatory expectations. Twenty-four months constitute both the industry benchmark and the threshold of reasonable evidence for regulatory justification without friction.

The spore itself could last much longer, but the integrated SCBI system cannot, and that distinction is what ultimately defines the number printed on the label.

---

## References

- ISO 11138-1:2017. Sterilization of health care products — Biological indicators — Part 1: Requirements for biological indicators.
- ISO 11138-3:2017. Sterilization of health care products — Biological indicators — Part 3: Biological indicators for moist heat sterilization processes.
- ICH Q1E: Evaluation of Stability Data. International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use.
- Russell, A.D. (2003). Bacterial spores and chemical sporicidal agents. *Clinical Microbiology Reviews*.
- Setlow, P. (2006). Spores of *Bacillus subtilis*: their resistance to and killing by radiation, heat and chemicals. *Journal of Applied Microbiology*.
- 3M Health Care (2021). Attest™ Rapid Readout Biological Indicator — Technical Manual.
- Mesa Labs (2022). SGM Biological Indicator System — Product Technical Documentation.

TerraGene

