

The Importance of Protein Quantification in Hygiene




Chemdye®
PRO1 MICRO

Abstract

Protein quantification is an essential tool for ensuring the effectiveness of cleaning processes in critical environments, such as hospitals and the food industry. Proteins, present in biological fluids and organic residues, can be a key indicator of contamination and hygiene failures if not properly removed. This document discusses the most commonly used methods for protein quantification and presents the Pro1 Micro system as an advanced example for real-time cleaning control and monitoring, validating the removal of proteins from surfaces.

1. Introduction

Cleaning and disinfection in critical areas, such as hospitals and food processing plants, require thorough control to prevent infections and ensure user safety. Protein residues on medical instruments or surfaces can act as a substrate for the proliferation of pathogenic microorganisms, compromising the effectiveness of sterilization and disinfection procedures. For this reason, protein quantification in hygiene processes has become a key standard to ensure effective cleaning.



2. Methods for Protein Quantification

Protein monitoring on surfaces after cleaning can be performed using several techniques:

- **Biuret Method:** A chemical method that detects peptide bonds in proteins, generating a color change that can be measured to quantify concentration.
- **Bradford Method:** This method uses a dye sensitive to proteins that changes color upon binding, allowing the measurement of small concentrations with high sensitivity.
- **Enzymatic Tests and Rapid Detection Kits:** These kits are designed for on-site testing, providing results in minutes. They are widely used for daily validation of cleaning processes in hospitals and laboratories.

Each method has its advantages and limitations, but the selection of the most suitable one depends on the environment and monitoring needs. In many cases, combinations of these methods are used to ensure surfaces are truly free of protein residues.



3. Pro1 Micro: An Advanced System for Cleaning Control

The **Pro1 Micro** is an innovative system that allows precise protein quantification to validate cleaning processes in critical environments. This device uses technology for rapid, quantitative detection of protein residues, providing real-time results.



Key features of the Pro1 Micro include:

- **Advanced Sensitivity:** Capable of detecting protein residues at extremely low concentrations, ensuring even the smallest traces are identified and removed during cleaning procedures.
- **Rapid Results:** Unlike other methods that may take hours, the Pro1 Micro delivers near-instantaneous data, allowing immediate adjustments to hygiene protocols.
- **Ease of Use:** Designed for busy environments, its swab design allows access to challenging areas during instrument washing.

The **Pro1 Micro** not only quantifies proteins but also helps users identify critical areas requiring additional attention, contributing to the continuous improvement of hygiene procedures.

4. Applications in Critical Environments

Protein quantification has direct applications in various critical environments, especially in hospitals and the food industry:

- **Hospitals and Clinics:** In the sterilization of medical instruments, the presence of proteins is a clear indicator of insufficient cleaning. The Pro1 Micro can verify the effectiveness of cleaning cycles before proceeding with sterilization, preventing the possibility of cross-contamination.
- **Food Industry:** In food processing plants, protein quantification ensures that equipment and surfaces are free of biological residues, reducing the risk of contamination from pathogens such as *E. coli* or *Listeria*.

5. Regulatory Standards and Best Practices

In critical environments such as hospitals, laboratories, and food processing plants, regulatory standards are essential for ensuring the proper implementation of cleaning procedures. International organizations have developed guidelines to ensure proper hygiene monitoring, and protein quantification has become one of the most recommended tools. Some of the most stringent standards include:

ISO 15883: This international standard defines the requirements for washer-disinfectors used in reprocessing medical instruments. ISO 15883 requires tests to ensure cleaning effectiveness, recommending protein tests as one of the most effective hygiene monitoring tools.

HTM 01-01: This technical guide from the UK's National Health Service (NHS) outlines detailed procedures for the decontamination and reprocessing of medical equipment. HTM 01-01 emphasizes the importance of cleaning tests, including protein tests, to ensure the complete removal of biological residues.

ANSI/AAMI ST79: Developed by the Association for the Advancement of Medical Instrumentation (AAMI), this standard is a comprehensive guide on the sterilization of medical devices in healthcare facilities. In this standard, protein tests are highlighted as a key tool for validating cleaning procedures and ensuring effective residue removal from surfaces and equipment.

These standards not only recommend protein quantification as part of hygiene monitoring but also promote its use as a primary tool for ensuring that cleaning protocols meet the necessary requirements to maintain safe and contamination-free environments.

The **Pro1 Micro** fits perfectly into these regulatory recommendations, providing a fast and accurate protein quantification system that ensures compliance with the highest international safety and hygiene control standards.

6. Conclusions

Protein quantification is an indispensable tool in hygiene control programs.

Its implementation ensures the effectiveness of cleaning processes and prevents risks associated with contamination. The **Pro1 Micro** is a standout example of advanced technology that not only quantifies proteins quickly and accurately but also offers continuous, real-time control, making it an ideal solution for environments where safety and hygiene are paramount.

In summary, protein quantification is not just a quality control measure but also ***an essential strategy for preventing infections and ensuring compliance*** with the highest safety standards in critical areas.

References

- Bradford, M. M. (1976). A rapid and sensitive method for the quantitation of microgram quantities of protein utilizing the principle of protein-dye binding. *Analytical Biochemistry*, 72(1-2), 248-254.
- Segal, J., & Muller, P. (2005). Protein quantification as an indicator of hygiene in medical equipment sterilization. *Journal of Hospital Infection*, 60(1), 15-20.
- World Health Organization (2016). Decontamination and reprocessing of medical devices for health-care facilities. *WHO Press*.
- International Organization for Standardization (ISO) (2018). ISO 15883-1:2018 Washer-disinfectors – General requirements, terms, and definitions and tests.
- U.S. Food and Drug Administration (2020). Guidelines for Validation of Cleaning Processes in Medical Device Manufacturing.
- European Committee for Standardization (CEN) (2017). EN 285: Large steam sterilizers
- Terragene. Pro1 Micro: Rapid Protein Quantification System.
- Terragene. Pro1 Micro Validation Report: Ensuring Cleanliness in Critical Environments.
- Smith, L., & Davis, P. (2021). Use of protein quantification in preventing hospital-acquired infections. *Healthcare Hygiene Journal*.
- Design and validation of a method for evaluating medical device cleanliness.