

A blue-tinted photograph of a group of medical professionals in white scrubs walking down a modern hospital hallway. In the foreground, a woman and a man are looking at a tablet together. Other staff members are visible in the background. The hallway has a glass railing on the left and a glass wall on the right.

Products & Cross Reference for Sterile Processing Department

Rev. 2 | April 2021



Index

Biological indicators

- 5** Bionova® BT96
- 6** Bionova® BT110
- 7** Bionova® BT220
- 8** Bionova® BT222
- 9** Bionova® BT223
- 10** Bionova® BT224

Process Challenge Devices | PCD

- 12** Bionova® KPCD220-C
- 13** Bionova® KPCD220-2
- 14** Bionova® KPCD222-C
- 15** Bionova® KPCD222-2
- 16** Bionova® KPCD224-C
- 17** Bionova® KPCD224-2
- 18** Bionova® PCD26-C

Bowie-Dick Test Packs

- 20** Chemdye® BD125X/1

Type 1 | Chemical indicators

- 22** Cintape® CT22
- 23** Cintape® CT40
- 24** Chemdye® CD42

Type 4 | Chemical indicators

- 26** Chemdye® CD16
- 27** Chemdye® CD29
- 28** Chemdye® CD40

Type 5 | Chemical indicators

- 30** Integron® IT12
- 31** Integron® IT26-C
- 32** Integron® IT26-1YS

Washing indicators

- 34** Chemdye® CDWU-Z
- 35** Chemdye® Splat CDWA4
- 36** Chemdye® CDWAH
- 37** Chemdye® CDWAH-U

Biological Indicators



Bionova® BT96

- Self-Contained Biological Indicator for VH2O2 sterilization monitoring | 30 minute readout.



Compatible with Bionova® Autoreaders
IC1020FR | IC1020FRLCD | MINIBIO



 **30' Super Rapid Readout**

VH2O2

Cross reference

Attest 1295 Celerity LCB044 Sterrad Velocity 43210




Bionova® BT110

- Self-Contained Biological Indicator for Ethylene Oxide sterilization monitoring | 4 hour readout.



Compatible with Bionova® Autoreaders
IC1020FR | IC1020FRLCD | MINIBIO



Cross reference

Attest 1294




Bionova® BT220

- Self-Contained Biological Indicator for Steam sterilization monitoring | 3 hour readout.
- Vacuum assisted and gravity air-displacement cycles at 121-135 °C.



Compatible with Bionova® Autoreaders
IC1020FR | IC1020FRLCD | MINIBIO



Cross reference

Attest 1292



Bionova® BT222

- Self-Contained Biological Indicator for Steam sterilization monitoring | 1 hour readout.
- Vacuum assisted and gravity displacement cycles 121-135 °C.



Compatible with Bionova® Autoreaders
IC1020FR | IC1020FRLCD | MINIBIO



1h Readout

STEAM

Cross reference

Verify Assert LCB031

Attest 1492V



Bionova® BT223

- Self-Contained Biological Indicator for IUSS monitoring (Flash Cycles) | 30 minutes readout.
- Gravity-air displacement cycles at 132-135 °C.



Compatible with Bionova® Autoreaders
IC1020FR | IC1020FRLCD | MINIBIO



Cross reference

Attest 1291

Attest 1492V



Bionova® BT224

- Self-Contained Biological Indicator for Steam sterilization monitoring | 20 minute readout.
- Vacuum assisted and gravity displacement cycles 132-135 °C.



Compatible with Bionova® Autoreaders
IC1020FR | IC1020FRLCD | MINIBIO



 **20'** Ultra Rapid
Readout

STEAM

Cross reference

Attest 1492V

Celerity LCB048



Process Challenge Devices | PCD



Bionova® KPCD220-C



- _25 PCDs + 25 SCBIs (positive control)
- _25 PCDs + 5 SCBIs
- _16 PCDs + 4 SCBIs

- Steam Test Pack PCD Kit Process Challenge Device for Steam sterilization processes | 3 hour readout.
- Vacuum assisted and gravity-air displacement between 121-135 °C.



Moving Front
Integrator



 **3h** Readout

STEAM

Cross reference

Attest 41382

Attest 41382F



Bionova® KPCD220-2



- _25 PCDs + 25 SCBIs (positive control)
- _25 PCDs + 5 SCBIs
- _16 PCDs + 4 SCBIs

- Steam Test Pack PCD Kit Process Challenge Device for Steam sterilization processes | 3 hour readout.
- Vacuum assisted and gravity-air displacement between 121-135 °C.



Unique Point Integrator



3h Readout

STEAM

Cross reference

Attest 41382

Attest 41382F



Bionova® KPCD222-C



- _25 PCDs + 25 SCBIs (positive control)
- _25 PCDs + 5 SCBIs
- _16 PCDs + 4 SCBIs

- Steam Test Pack PCD Kit Process Challenge Device for Steam sterilization processes | 1 hour readout.
- Vacuum assisted and gravity-air displacement between 121-135 °C.



Moving Front
Integrator



 **1h** Readout

STEAM

Cross reference

Verify Assert LCB032 | LCB033 | LCB043

Attest 1496V | 1492VF



Bionova® KPCD222-2



_25 PCDs + 25 SCBIs (positive control)
_25 PCDs + 5 SCBIs
_16 PCDs + 4 SCBIs

- Steam Test Pack PCD Kit Process Challenge Device for Steam sterilization processes | 1 hour readout.
- Vacuum assisted and gravity-air displacement between 121-135 °C.



Cross reference

Verify Assert LCB032 | LCB033 | LCB043

Attest 1496V | 1492VF



Bionova® KPCD224-C



- _25 PCDs + 25 SCBLs (positive control)
- _25 PCDs + 5 SCBLs
- _16 PCDs + 4 SCBLs

- Steam Test Pack PCD Kit Process Challenge Device for Steam sterilization processes | 20 minute readout.
- Vacuum assisted and gravity-air displacement between 132-135 °C.



Moving Front
Integrator



20'Rapid
Readout

STEAM

Cross reference

Celerity LCB049 | LCB050
LCB053 | LCB033

Attest 1496V

Attest 41482V
41482VF



Bionova® KPCD224-2



- _25 PCDs + 25 SCBLs (positive control)
- _25 PCDs + 5 SCBLs
- _16 PCDs + 4 SCBLs

- Steam Test Pack PCD Kit Process Challenge Device for Steam sterilization processes | 20 minute readout.
- Vacuum assisted and gravity-air displacement between 132-135 °C.



Cross reference

Celerity LCB049 | LCB050
LCB053 | LCB033



Attest 1496V



Attest 41482V
41482VF



Bionova® PCD26-C

- Process Challenge Device Test pack for Steam sterilization processes with moving front integrator.
- Monitoring dynamic air removal and gravity displacement between 121-135 °C.



Cross reference

Attest 41360



Bowie-Dick Test Packs



Chemdye® BD125X/1

- Bowie-Dick Test Pack Chemical Indicator for Steam sterilization. With warning sheet.
- Testing air removal efficiency for cycles of 132 °C for 4 min. and 134 °C for 3.5 min.
- Turns from Yellow to Dark Brown/Black. 100% Lead Free.



Cross reference

Comply 00135LF

Verify EQC009




Type 1 Chemical Indicators



Cintape® CT22

- Steam Indicator Tape | Type 1 Chemical Indicator | 18 mm.
- Vacuum assisted and gravity-air displacement between 121-135 °C.
- Initial color:  | Final color:  



Cross reference

Comply 1322-18mm

Browne 0142 | 0143 | 04147

Verify 260010



Cintape® CT40

- VH2O2 Indicator Tape | Type 1 Chemical Indicator | 20 mm.
- Initial color:  | Final color: 



Cross reference

Comply 1228

Browne 0142 | 0143 | 04147

Verify 260010



Chemdye® CD42

- Process Indicator for monitoring Plasma or Vaporized Hydrogen Peroxide sterilization processes.



Cross reference

Comply 1248



Celerity PCC074



Verify PCC061



ASP 14100



Type 4 Chemical Indicators



Chemdye® CD16


- Multivariable Process Indicator for monitoring Ethylene Oxide sterilization processes.
- Conditions: 600 mg/l EO, RH 60%, 30 °C, 30 min. | 600 mg/l EO, RH 60%, 54 °C, 20 min.



Cross reference


Comply 1251

Browne 2563



Chemdye® CD29

- Multivariable Process Indicator for monitoring Steam sterilization processes.
- Conditions: 121 °C, 15 minutes | 134 °C, 3.5 minutes.



Cross reference

Cross-checks CI-111 | 107 | 134 | 139

Comply 1250

Browne 2551



Chemdye® CD40

- Multivariable Process Indicator for monitoring Plasma or Vaporized Hydrogen Peroxide sterilization processes.
- Conditions: 6 minutes, 50 °C, 2,3 mg/l H₂O₂.



Cross reference

CVP2-214-241, 214-243
CVP4 214-231, 214-233
CVP6 214-221, 214-223

Browne 2563

Verify PCC048



Type 5 Chemical Indicators



Integron® IT12

- Two level Chemical Integrator for Ethylene Oxide sterilization processes.




Cross reference

Browne 2420



Integron® IT26-C

- Type 5 Integrating Indicator. Size: 52 x 21 mm.




Cross reference


Comply Sterigage 1243A | 1243B | 1243RE



Comply Thermalog 2134 mm




Verify PCC064 | PCC064B
PCC065 | PCC065B | PCC067
PCC067B | PCC039



Integron® IT26-1YS

- Unique point Chemical Integrator for Steam sterilization processes.



Cross reference

Verify 800510 long | 801010 short

CSP5-SV1 211-224 | 211-225 | 211-226



Washing Indicators



Chemdye® CDWU-Z

- Ultrasonic Cleaning Indicator for ultrasonic cavitation performance tests.



Cross reference

Sonno Check TI108



Chemdye® Splat CDWA4

- Cleaning Indicator for cleaning performance tests.



Cross reference

Verify LCC010 | LCC012

ProFormance Washer Test TOSI



Chemdye® CDWAH

- Holder for Chemdye® Splat for Monitoring Washer Disinfectors.



Cross reference

Verify LCC010



Chemdye® CDWAH-U

- Holder for Chemdye® Splat for Monitoring Ultrasonic Washing Machines.





Let's work together
to create a better future.

