

Artificial intelligence helping you make better decisions.

Terragene® has developed a simple and accurate analysis device for chemical and washing indicators. Get clear results with precision that goes beyond human capabilities. The scanner works together with proprietary Bionova® Cloud software allowing to integrate all your results in the same platform: washing indicators, chemical indicators, protein-based hygiene monitoring system and biological indicators.



Smart solutions

AI technology assists in making more accurate interpretations.



Digitize all your results

All your indicators together in the same Bionova® Cloud platform: the best choice quality management inside the SPD/CSSD. No more physical databases.



Whole traceability

All the information is automatically linked and generated: indicator SKU, indicator lot number, expiration date, final result and sterilizer/washer where it was exposed.



Performance you can rely on

Outstanding hardware and software technology working together to give accurate information avoiding any human potential mistake from Outstanding hardware and software technology working together to give accurate information avoiding any human potential mistake going from a missed indicator, to the result interpretation and even with the equipment monitoring.



Result imaging is the first step in building a digital transformation plan.



Getting easier and faster access to reliable information across multiple facilities

Trazanto now allows the automatic incorporation of results of chemical indicators and washing indicators in the same traceability system where the rest of the Terragene® monitoring indicators are automatically managed and analyzed. Terragene® devices allows to manage all together in the same platform: Bionova® biological indicators, Chemdye® ProMicro (hygiene monitoring system based on the absolute quantification of total proteins), Chemdye® chemical indicators, Chemdye® washing Integron® chemical indicators.

- LARGE SURGERY CENTERS
- HOSPITALS
- SMALL SURGERY CENTERS
- COMMUNITY HEALTH SERVICES
- CLINICS
- DENTAL OFFICES

- BIOLOGICAL INDICATORS READOUT
- HYGIENE MONITORING PEN READOUT
- CHEMICAL AND WASHING INDICATORS READOUT


Intelligent algorithms developed for automated reading of Terragene® indicators

Automatic system for quality control and traceability of cleaning, disinfection and sterilization processes

Fully compatible with **Bionova® Cloud traceability software.**

Usage
For readout and detection of chemical indicators

COMPATIBLE PRODUCTS



Two main modules:

- Hardware module composed by the Trazanto optical reader.
- Bionova® Cloud software-implemented analysis module.



Disruptive technology for sterilization process management

